

# AirCard™ Wireless Phone Capability

Now an option in ECO™-14, Vantage™ Plus, Vanguard™ Plus and Pathfinder™ Plus



*AirCard® wireless phone capability enables your PRC communication aid to act as a mobile phone. The AirCard inserts into the PCMCIA slot to allow the client to send and receive voice calls and text messages.*

## Capabilities

*Send and receive phone calls entirely through your PRC communication aid.*

ECO-14, Vantage Plus, Vanguard Plus and Pathfinder Plus from PRC are all now compatible with the Sierra AirCard® technology.

The AirCard utilizes the SIM card from a standard GSM cellular phone and simply plugs into the PCMCIA slot in the PRC device.

This eliminates the need for a cable connection between your phone and your communication aid. And you no longer need to charge both your device and your cell phone. The AirCard is powered through the device.

You simply select a cell phone plan from a supporting GSM network in your area. For a list of supporting networks and compatible AirCards, contact the PRC Service Department at 1-800-262-1990.



The AUDIO CABLE adapts the AirCard to your PRC device.

Antenna provided with the AirCard.

*Send and receive text messages entirely through your PRC communication aid.*

For text messaging you can formulate your messages using Unity®, WordCore™ or any other vocabulary program compatible with our ECO-14, Vantage, Vanguard and Pathfinder devices.

The mobile phone function allows you to send and receive calls, listen to the caller and talk entirely through your PRC communication aid.

# AirCard Wireless Phone Capability

## What do I need to make and receive calls and for text messages?

To make and receive calls and use text messaging you need

- the AirCard with antenna,
- a cell phone service contract for a GSM network (including a SIM\* for the AirCard)
- an AUDIO CABLE to make and receive calls.

*PRC devices released after October 1, 2005 have all the necessary hardware updates for phone capability. Devices manufactured prior to October, 2005 will require hardware updates and must come to PRC headquarters to be updated. Call PRC Sales office at (800) 262-1933 for more information on product upgrades and pricing.*

### Simple Steps

- 1) Get your PRC ECO-14, Pathfinder, Vanguard or Vantage.
- 2) Buy an AirCard accessory package from PRC.
- 3) Get a GSM phone plan.\*\*
- 4) Insert the SIM into the AirCard.
- 5) Insert the AirCard into your PRC device.
- 6) START TALKING!

*ECO-14 AirCard Accessory (complete, includes card)*

AC-E14 Price: \$495

*Pathfinder AirCard Accessory (complete, includes card)*

AC-PF Price: \$495

*Vanguard AirCard Accessory (complete, includes card)*

AC-VG Price: \$495

*Vantage AirCard Accessory (complete, includes card)*

AC-VT Price: \$495

The AirCard features allow for cell phone communications, which is a must for any communication device user. One of my favorite features is having the ability to send text messages, which is sometimes useful in louder places like sporting facilities.



I feel safer and much more independent knowing I have a cell phone I can access.

*Anthony Arnold*

\* A SIM (Subscriber Identification Module) is provided through a wireless service provider. The SIM contains the cell phone number and other network activation information.

\*\*PRC does not provide wireless service and customers will have to provide their own mobile phone service account. The AirCard currently works with any service provider that supports the GSM network (e.g., AT&T, T-Mobile). Customers will be responsible for their own phone charges and PRC cannot accept responsibility for problems of mobile phone service provision.

**To order a ECO-14, Pathfinder, Vanguard or Vantage or to get a price quote, call toll-free 1-800-262-1933 or visit our online store at [www.prentrom.com](http://www.prentrom.com)**

ECO-14, Pathfinder, Vanguard, Vantage, and WordCore are trademarks of Prentke Romich Company · Unity® is a registered trademark, in the USA, UK, and Canada, and the property of Semantic Compaction Systems · AirCard is a registered trademark of Sierra Wireless · AT&T is a registered trademark of AT&T, Inc. · T-Mobile is a registered trademark of Deutsche Telekom AG.



PRC · 1022 Heyl Road · Wooster, Ohio 44691  
ACF Rev. 08-13-08