

ECO™ -14



Standard Dedicated Device Configuration

(E14-AEN-STD) \$7,595

Dedicated device meets Medicare/Medicaid funding requirements.

Standard Configuration includes:

- Large (14.1") XGA TFT display (1024 x 768)
- 8 to 9 hours of battery life
- Built-in DVD/CD player and CD burner
- Integrated Bluetooth® connectivity (serial port protocol only)
- External monitor port (XGA)
- Mounting plate
- Two-position table-top stand
- 512MB RAM, 500Mhz, 60GB storage
- USB keyboard & mouse
- Three USB2 ports
- Integrated speaker
- Integrated keyboard
- Access methods: Touch screen, on-screen keyboard, keyboard and mouse, visual and auditory scanning, 1-, 2- or multi-switch scanning, wheelchair joystick, headmouse or pointer
- Environmental controls
- Built-in infrared capabilities
- Operator manuals on CD
- PRC service and support including a one-year warranty
- All versions of the word-based WordCore™ vocabulary system
- Symbols/language options: Unity 4, 8, 15, 45, 60, 84, 144, Adult Acquired Module
- Language representation methods supported: Single meaning pictures, spelling, word prediction, orthographic word selection, Semantic Compaction
- Language Activity Monitor
- Customizable block scanning
- PRC Application and Support Software (PASS™) with added features

Integrated Device Configuration

(E14-AEN-INTG) \$7,595

By selecting this Integrated version, the system allows you to switch between communication and computer features.

Integrated Configuration includes:

- All Standard Configuration features shown above
- "One touch" transition from software programs to speech
- Windows XP® Home Edition
- Use Unity® as your onscreen keyboard for Windows XP programs
- CD RW/DVD-ROM combo drive

Color Options

Select one of the following color options:

Ocean Blue (BLU), Midnight Black (BLK), Coral Pink (PNK) or Sunset Orange (ORG)

A black case will be shipped if no color choice is made at the time of the order. A service charge will be applied on units returned to PRC for a color change.



PRC · 1022 Heyl Road · Wooster, OH 44691
Main Phone: (800) 262-1984 · Email: info@prentrom.com
Sales Phone: (800) 262-1933 · Web: www.prentrom.com

Other Product Options

These options must be specified at the time of purchase.

Integrated Headpointing (Tracker™ technology)
(HP) add \$895

IntelliTools® Classroom Suite Bundle
(ICS) add \$299

PCS for ECO-14
(PCS) add \$225
5,000 PCS symbols from Mayer-Johnson pre-loaded into ECO-14 for your use.

WordPower™ Version - Available Late Spring 2007
(E14-AEN-WDPR) add \$350



Upgrading your ECO-14

(PCS) (after initial sale) add \$250
For more information on available upgrades, visit our online store at www.prentrom.com or contact our Sales Department at (800) 262-1933.

Accessories

144-location Keyguard Kit (144KG-E14)..... \$ 135
84-location Keyguard Kit (84KG-E14) \$ 135
60-location Keyguard Kit (60KG-E14) \$ 135

AirCard™ (complete with cable and card) (AC-E14) \$ 495
Battery Charger (BC-E14) \$ 75
Carrying Case (CC-E14) \$ 125
Daessy™ Desktop Mount (DDM-1) \$ 410



Extended Warranty*

	Order Number	One-Year Package	Two-Year Package	Three-Year Package	Four-Year Package	Five-Year Package
Gold	EW-EC14-G	\$465	\$888	\$1269	\$1584	\$1890
Platinum	EW-EC14	\$600	\$1164	\$1665	\$2100	\$2487

*Excludes USB Keyboard & Mouse.

When not specified, pricing includes coverage at the Platinum level of service. A five-year package may be purchased only before the first year of warranty expires. Extended Security may be purchased in increments of one to five years. If the warranty has expired, the Extended Security may be re-instituted if the device is inspected by service and repaired at the expense of the customer. After the device is re-certified, additional years of Extended Security may be purchased to cover the device up to six years from the purchase date. Please include applicable state sales tax in the following states, or include a tax exemption certificate: AR, CT, HI, IA, KS, LA, MS, SC, SD, TX, UT and WV. Warranty not available for sale in the state of Virginia.

ECO-14 Specifications

Size 13.2" x 11.2" x 2.9"
Weight 6 lbs. 13 oz. (with single, swappable battery pack)
Battery Lithium Ion, 8-9 hours of battery life
LCD Screen Size 14.1"
LCD Screen Resolution 1024 x 768
Touchscreen Sensitivity Finger & stylus
Synthesized Speech RealSpeak, DECtalk and AT&T
Digitized Speech Record Time: 660K/minute



ECO, Language Activity Monitor (LAM), PASS, Pathfinder, SpringBoard, Vanguard, Vantage, and WordCore are trademarks of Prentke Romich Company · Unity® is a registered trademark, in the USA, UK, and Canada, and the property of Semantic Compaction Systems · Semantic Compaction is a trademark of Semantic Compaction Systems · AirCard is a registered trademark of Sierra Wireless · AT&T Natural Voices is a trademark of AT&T Corporation · Bluetooth Wireless Technology is a registered trademark of Bluetooth SIG, Inc. · Daessy is a trademark of Daedalus Technologies, Inc. · Fonix-DECtalk Text-to-Speech software is a copyright of Fonix Corporation · IntelliTools is a registered trademark of IntelliTools, Inc. · RealSpeak is a product of Scansoft · Tracker is a trademark of Madentec · Windows and Windows XP are registered trademarks of Microsoft Corporation · WordPower is a trademark of Inman Innovations, Inc. · WordQ is a trademark of Bloorview MacMillan Children's Centre.
Product prices and features are subject to change without notice.