





Comparison of SpringBoard Plus and SpringBoard Lite

	SpringBoard Plus 	SpringBoard Lite 
Vocabularies	8, 15, and 32 location word base vocabulary in child and adult male and female voices. 4 and 8 location phrase base vocabulary in child and adult male and female voices. There is no difference between vocabulary recorded in the child's voice or adult male and female voices.	32, 15, and 8 location word-based (child) 8 and 4 location phrase based (child) 8, 15, and 32 location word-based adult vocabulary (male and female) 8 and 4 location phrase-based vocabulary (male and female); 15 Word Cues and CVA
Key size	32 location 3/4" X 3/4" 15 location 1 3/16" X 1 1/8" 8 location 1 1/2" X 1 3/4" 4 location 3" X 1 3/4" All Activity rows 3/4" X 3/4"	32 location 11/16" X 3/4" 15 location 1 3/16" X 7/8" 8 location 1 2/3" X 1 3/8" 4 location 3" X 1 3/8" All Activity rows 11/16" X 3/4"
Screen Size	6-1/8" x 4-5/8"	6-1/16" x 3- 5/8"
Weight	3 lbs 7oz	2 lbs 8 oz
Access	Direct select Mouse 1 and/or 2 switch IR Headpointing (optional)	Direct select Mouse 1 and/or 2 switch
Keyguards	Set of 8, 15, and 32 (clear)	Set of 8, 15, and 32 8 and 15 are black 32 clear
Voice Options	prerecorded child, male and female	prerecorded child, male and female
Speech capabilities	All recorded speech; has the capability for additional 200 minutes.	All recorded speech; has the capability for additional 500 minutes.
Graphic handling options	Converts all types of graphics to 256 color (BMP 256 quality)	Converts all types of graphics to a 18K color (Jpeg like quality)

	SpringBoard Plus	SpringBoard Lite
		
Visual Scenes	no	yes
LCD	Screen is harder to see at angles	Brighter display – can be seen at different angles
Stand	Standard black plastic stand; Wire stand additional	Standard wire stand
Carrying Case	Optional Carrying case with plastic covering the screen.	Optional Fitted foam carrying case with removable cover.
USB Jumpdrive (for backing up the device)	Standard	Jump drive standard
USB cable	Standard	Optional USB Cable (regular USB cable).
Mounting Plate	Standard	Optional
Bluetooth for switches and computer access	no	Optional \$200 for Bluetooth inside; \$250 for the receiver Bluetooth to computer.
Colors	red	Choice of seafoam green, black, white, sunset orange, or vineyard grape
Toolbox	Full option only	Simple and full toolbox standard
Languages available	English Bilingual (Spanish/English)	English Bilingual (Spanish/English)
Ability to program on computer	yes via PASS for SpringBoard Plus	yes via PASS for SpringBoard Lite
PCS icons	\$395 after purchase	Optional \$375 at time of purchase \$395 after purchase
Cost	No longer available.	\$2495 English \$2745 Bilingual

Terese Phelps M.Ed., PRC Regional Consultant 02/16/09 SB-SBL