



Download Directions

Prentke Romich Application and Support Software (PASS)

1. GO TO www.prentrom.com



Find Your Consultant Support » PASS Demo Software

ZIP


PASS Demo Software

Current Devices

- Essence
- Vantage Lite
- SpringBoard Lite
- Vanguard Plus
- ECO2

Previous Devices

- ECO-14



We've gone digital!
Look for our e-book selection at Amazon.com and BarnesandNoble.com

[Read more . . .](#)

4. CHOOSE STUDENT'S DEVICE

Support » PASS Demo Software


PASS Demo Software


Current Devices

- Essence
- Vantage Lite
 - [LATEST \(Version 5.08\) Vantage Lite PASS - Supports device system software 5.08, 5.07, 5.06, or 5.05](#)
 - [LATEST \(Version 5.08\) Vantage Lite Spanish Bilingual PASS - Supports device system software 5.08, 5.07, or 5.05](#)
 - [\(Version 5.04.2\) Vantage Lite PASS Software Installation](#)
- SpringBoard Lite

5. CHOOSE PASS SOFTWARE INSTALLATION

Make sure software version matches one in device (look under Maintenance Menu in Toolbox on device).

 **File PASS Manual**
Size 2 MB

 **File: VTVG_PASS_AEN_INSTALL_5.08.exe**
Size: 272 MB

6. CLICK DOWNLOAD NOW

7. Follow the directions on your computer screen until download is complete.